

ENSAYOS DE LABORATORIO DE LA INDUSTRIA ALIMENTICIA

DE USO GENERAL GENERAL USE MEDIA	
Aeromonas Aeromonas	Staphylococcus Staphylococcus
<u>Base de Medio CromoCen AGN</u> <u>CromoCen AGN Medium Base</u>	<u>Agar Manitol Salado</u> <u>Manitol Salt Agar</u>
	<u>Base de Agar Baird-Parker</u> <u>Baird-Parker Agar Base</u>
Anaerobios Anaerobes	<u>Base de Agar Sangre</u> <u>Blood Agar Base</u>
	<u>m-Medio de Chapman</u> <u>m-Chapman Medium</u>
<u>Base de Agar Sangre</u> <u>Blood Agar Base</u>	
<u>Medio Leche de Crossley</u> <u>Crossley Milk Medium</u>	Enterococcus faecalis Enterococcus faecalis
<u>Medio Tioglicolato</u> <u>Thioglycollate Medium</u>	
	<u>Agar de MacConkey sin Cristal Violeta</u> <u>ni Sal</u> <u>MacConkey Agar without Crystal Violet</u> <u>and Salt</u>
Bacillus cereus Bacillus cereus	
	ANALISIS DE ALIMENTOS CONGELADOS FROZEN FOODS ANALYSIS
<u>Agar Dextrosa y Triptona</u> <u>Dextrose Tryptone Agar</u>	
<u>Agar Nutriente</u> <u>Nutrient Agar</u>	<u>Agar para Conteo en Placas</u> <u>Plate Count Agar</u>
<u>m-Caldo Glucosa y Triptona</u> <u>m-Glucose Tryptone Broth</u>	<u>Agar Dextrosa y Triptona</u> <u>Dextrose Tryptone Agar</u>
	<u>Agar Malta</u> <u>Malt Agar</u>
Campylobacter Campylobacter	<u>Agar Nutriente</u> <u>Nutrient Agar</u>
	<u>Agar Papa y Dextrosa</u> <u>Potato Dextrose Agar</u>

<u>Base de Agar Sangre Columbia</u> <u>Columbia Blood Agar Base</u>	<u>Agar Triptona Soya</u> <u>Tryptone Soy Agar</u>
<u>Base de Medio Brucella</u> <u>Brucella Medium Base</u>	<u>Agar Violeta Rojo Bilis</u> <u>Violet Red Bile Agar</u>
<u>Medio Tioglicolato</u> <u>Thioglycollate Medium</u>	<u>Base de Agar Endo</u> <u>Endo Agar Base</u>
	<u>Base de Medio CromoCen AGN</u> <u>CromoCen AGN Medium Base</u>
E. coli 0157:H7 E. coli 0157:H7	<u>Base de Medio CromoCen SC</u> <u>CromoCen SC Medium Base</u>
	<u>m-Base de Medio Endo</u> <u>m-Endo Medium Base</u>
<u>Agar de MacConkey con Sorbitol</u> <u>MacConkey Sorbitol Agar</u>	<u>Caldo Bilis Verde Brillante</u> <u>Brilliant Green Bile Broth</u>
<u>CromoCen ECCS</u> <u>CromoCen ECCS</u>	<u>Caldo de EC</u> <u>EC Broth</u>
	<u>m-Caldo Triptona Glucosa Extracto L</u> <u>m-Tryptone Glucose Extract Broth L</u>
Grupo coli-aerogenes Coli-aerogenes Group	<u>CromoCen CC</u> <u>CromoCen CC</u>
	<u>CromoCen ECCS</u> <u>CromoCen ECCS</u>
<u>Agar Citrato de Simmons</u> <u>Simmons Citrate Agar</u>	<u>m-Medio de Sabouraud</u> <u>m-Sabouraud Medium</u>
<u>Agar Desoxicolato</u> <u>Desoxycholate Agar</u>	<u>Medio Tioglicolato</u> <u>Thioglycollate Medium</u>
<u>Agar Eosina Azul de Metileno (Mod.)</u> <u>Eosin Methylene Blue Agar (Mod.)</u>	
<u>Agar de MacConkey</u> <u>MacConkey Agar</u>	ANALISIS DE CONSERVAS CANNED FOODS ANALYSIS
<u>Agar Tres Azúcares y Hierro</u> <u>Triple Sugar Iron Agar</u>	
<u>Agar Violeta Rojo Bilis</u> <u>Violet Red Bile Agar</u>	<u>Agar Dextrosa y Triptona</u> <u>Dextrose Tryptone Agar</u>
<u>Agua Triptona</u> <u>Tryptone Water</u>	<u>Agar Extracto de Levadura</u> <u>Yeast Extract Agar</u>
<u>Base de Agar Endo</u> <u>Endo Agar Base</u>	<u>Agar Extracto de Malta</u> <u>Malt Extract Agar</u>
<u>Base de Medio CromoCen SC</u> <u>CromoCen SC Medium Base</u>	<u>Agar de MacConkey</u> <u>MacConkey Agar</u>
<u>m-Base de Medio Endo</u> <u>m-Endo Medium Base</u>	<u>Agar Malta</u> <u>Malt Agar</u>
<u>Caldo Bilis Verde Brillante</u> <u>Brilliant Green Bile Broth</u>	<u>Agar Nutriente</u> <u>Nutrient Agar</u>

<u>Caldo de EC</u> <u>EC Broth</u>	<u>Agar Papa y Dextrosa</u> <u>Poptato Dextrose Agar</u>
<u>Caldo Lactosado</u> <u>Lactose Broth</u>	<u>Agar Triptona Soya</u> <u>Tryptone Soy Agar</u>
<u>Caldo Lauril Triptosa</u> <u>Lauril Tryptose Broth</u>	<u>Agar Violeta Rojo Bilis</u> <u>Violet Red Bile Agar</u>
<u>Caldo Púrpura de MacConkey</u> <u>MaConkey Broth (purple)</u>	<u>Base de Medio CromoCen SC</u> <u>CromoCen SC Medium Base</u>
<u>CromoCen CC</u> <u>CC CromoCen</u>	<u>m-Caldo Glucosa y Triptona</u> <u>m-Glucose Tryptone Broth</u>
	<u>m-Caldo Triptona Glucosa Extracto L</u> <u>m-Tryptone Glucose Extract Broth L</u>
Mohos y Levaduras Yeasts and Moulds	<u>CromoCen CC</u> <u>CromoCen CC</u>
	<u>Medio Leche de Crossley</u> <u>Crossley Milk Medium</u>
<u>Agar Dextrosa de Sabouraud</u> <u>Sabouraud Dextrose Agar</u>	<u>m-Medio de Sabouraud</u> <u>m-Sabouraud Medium</u>
<u>Agar Extracto de Malta</u> <u>Malt Extract Agar</u>	
<u>Agar Malta</u> <u>Malt Agar</u>	ANALISIS EN LA INDUSTRIA CERVECERA BEVERAGE INDUSTRY ANALYSIS
<u>Agar Maltosa de Sabouraud</u> <u>Sabouraud Maltose Agar</u>	
<u>Agar Papa y Dextrosa</u> <u>Potato Dextrose Agar</u>	<u>Agar Extracto de Levadura</u> <u>Yeast Extract Agar</u>
<u>Caldo Extracto de Malta</u> <u>Malt Extract Broth</u>	<u>Base de Medio CromoCen SC</u> <u>CromoCen SC Medium Base</u>
<u>m-Medio de Sabouraud</u> <u>m-Sabouraud Medium</u>	<u>Caldo Lactosado</u> <u>Lactose Broth</u>
	<u>Caldo Púrpura de MacConkey</u> <u>MaConkey Broth (purple)</u>
Organismos viables Viable organisms	<u>CromoCen CC</u> <u>CromoCen CC</u>
<u>Agar para Conteo en Placas</u> <u>Plate Count Agar</u>	ANALISIS DE PRODUCTOS LACTEOS DAIRY PRODUCTS ANALYSIS
<u>Agar Nutriente</u> <u>Nutrient Agar</u>	
<u>Agar Triptona Glucosa Extracto</u> <u>Tryptone Glucose Extract Agar</u>	<u>Agar Citrato de Simmons</u> <u>Simmons Citrate Agar</u>
<u>Agar Triptona Soya</u>	<u>Agar para Conteo en Placas</u>

<u>Tryptone Soy Agar</u>	<u>Plate Count Agar</u>
<u>Base de Agar Sangre</u> <u>Blood Agar Base</u>	<u>Agar Desoxicolato</u> <u>Desoxycholate Agar</u>
<u>m-Caldo Triptona Glucosa Extracto C</u> <u>m-Tryptone Glucose Extract Broth C</u>	<u>Agar Extracto de Levadura</u> <u>Yeast Extract Agar</u>
<u>m-Caldo Triptona Glucosa Extracto L</u> <u>m-Tryptone Glucose Extract Broth L</u>	<u>Agar Leche</u> <u>Milk Agar</u>
<u>Caldo Triptona Soya</u> <u>Tryptone Soy Broth</u>	<u>Agar Nutriente</u> <u>Nutrient Agar</u>
	<u>Agar Papa y Dextrosa</u> <u>Poptato Dextrose Agar</u>
Otras enterobacterias (No se los del grupo coli-aerogenes, Salmonella y Shigella) Other enteric bacteria (Coli-aerogenes group, Salmonella and Shigella not included)	<u>Agar Triptona Glucosa Extracto</u> <u>Tryptone Glucose Extract Agar</u>
	<u>Agar Triptona Soya</u> <u>Tryptone Soy Agar</u>
	<u>Agar Violeta Rojo Bilis</u> <u>Violet Red Bile Agar</u>
	<u>Base de Agar Endo</u> <u>Endo Agar Base</u>
<u>Agar Hierro de Kligler</u> <u>Kligler Iron Agar</u>	<u>Base de Agar Sangre Columbia</u> <u>Columbia Blood Agar Base</u>
<u>Agar de MacConkey sin Cristal Violeta ni Sal</u> <u>MacConkey Agar withoutCrystal Violet and Salt</u>	<u>Base de Medio Brucella</u> <u>Brucella Medium Base</u>
<u>Agar Tres Azúcares y Hierro</u> <u>Triple Sugar Iron Agar</u>	<u>Base de Medio CromoCen SC</u> <u>CromoCen SC Medium Base</u>
<u>Base de Agar Selectivo para Yersinia</u> <u>Yersinia Selective Agar Base</u>	<u>m-Base de Medio Endo</u> <u>m-Endo Medium Base</u>
<u>Base de Agar Urea</u> <u>Urea Agar Base</u>	<u>Caldo Bilis Verde Brillante</u> <u>Brilliant Green Bile Broth</u>
<u>Caldo Triptona Soya</u> <u>Tryptone Soy Broth</u>	<u>Caldo Brucella</u> <u>Brucella Broth</u>
	<u>Caldo Corazón</u> <u>Heart Infusion Broth</u>
Pseudomonas Pseudomonas	<u>Caldo de EC</u> <u>EC Broth</u>
	<u>Caldo Lactosado</u> <u>Lactose Broth</u>
<u>Base de Medio CromoCen AGN</u> <u>CromoCen AGN Medium Base</u>	<u>Caldo Lauril Triptosa</u> <u>Lauril Tryptose Broth</u>

	<u>m-Caldo Triptona Glucosa Extracto L</u> <u>m-Tryptone Glucose Extract Broth L</u>
Salmonella y Shigella <i>Salmonella and Shigella</i>	<u>Caldo Triptosa Fosfato</u> <u>Tryphtose Phosphate Broth</u>
	<u>CromoCen CC</u> <u>CromoCen CC</u>
<u>Agar Desoxicolato y Citrato</u> <u>Desoxycholate Citrate Agar</u>	<u>Medio Leche de Crossley</u> <u>Crossley Milk Medium</u>
<u>Agar Entérico de Hektoen</u> <u>Hektoen Enteric Agar</u>	<u>m-Medio de Sabouraud</u> <u>m-Sabouraud Medium</u>
<u>Agar Hierro de Kligler</u> <u>Kligler Iron Agar</u>	<u>Base de Medio CromoCen AGN</u> <u>CromoCen AGN Medium Base</u>
<u>Agar Lisina y Hierro</u> <u>Lysine Iron Agar</u>	<u>CromoCen ECCS</u> <u>CromoCen ECCS</u>
<u>Agar S.S. (Mod.)</u> <u>S.S. Agar (Mod.)</u>	
<u>Agar Sulfito Bismuto</u> <u>Bismuth Sulphite Agar</u>	
<u>Agar Tres Azúcares y Hierro</u> <u>Triple Sugar Iron Agar</u>	
<u>Agar Verde Brillante</u> <u>Brilliant Green Agar</u>	
<u>Agar X.L.D.</u> <u>X.L.D. Agar</u>	
<u>Agua de Peptona Buferada</u> <u>Buffered Peptone Water</u>	
<u>Base de Caldo Selenito</u> <u>Selenite Broth Base</u>	
<u>Base de Caldo Tetrionato</u> <u>Tetrathionate Broth Base</u>	
<u>Base de Medio CromoCen SC</u> <u>CromoCen SC Medium Base</u>	
<u>Base de Medio para Descarboxilación</u> <u>de Aminoácidos</u> <u>Aminoacids Descarboxylation Medium</u> <u>Base</u>	
<u>Caldo GN</u> <u>GN Broth</u>	

Fuente: <http://www.biocen.cu/producto/indicemc/alimen.htm>